

Classification of NYC Aerosols by X-Ray and Optical Methods

By: John Sangobowale (Mount St. Michael Academy)
and William Dennis (John F. Kennedy H.S.)

Mentors: Marc Cesaire (Graduate student) and
Dr. Elizabeth Rudolph

EAS Department CCNY

Purpose

The Principal objective of our work is:

- Elemental Characterization of Aerosols collected by two methods for comparison:
 - EBAM: beta mass attenuation
 - Millipore apparatus: Vacuum Filtration
- Ultimately: to understand how weather patterns affect the chemical composition and darkness of aerosol particles

What are Aerosols?

- Aerosols are small solid or liquid particles suspended in the atmosphere. Their sizes vary from a few nanometers (0.000000001 meters) to almost 100 micrometers (0.0001 m, the thickness of a hair.
 - Volcanic dust
 - Combustion products
 - Soot
 - Smoke

Origin of Aerosols

- Aerosols originate both from natural and man-made (anthropogenic) sources . They can be directly emitted as particles (primary aerosols) or they can also be the result of chemical reactions (secondary aerosols).


Why do we study aerosols?

- Health issues: direct relationship between lung diseases


- Asthma
- Lung cancer

and particles 2.5 microns and smaller

- Safety and Security: Chemical classification of air masses to be able to recognize unusual or dangerous new patterns




Collection: Environmental Beta Attenuation Mass Monitor (EBAM)



2.5 micron particles

Analytical Methods

- X-RAY Fluorescence Spectrometry (XRF)



Philips PW 1400 X-Ray Fluorescence Spectrometer (XRF)





Classification by Element Constituents

- Using chemistry and elemental analysis: In particular the following elements:
 - Titanium (Ti)
 - Iron (Fe)
 - Sulfur (S)
- Logic for element selection:
 - Preliminary scans across spectra for all elements shows titanium iron and sulfur in intensities much above background
 - Sulfur compounds are produced through burning of fuels
 - Many elements are potentially harmful to humans especially in particles of 2.5 microns or less


XRF Results



XRF Results June 2004



XRF Results July 2004



Characterization by optical properties


- Correlating sample optical properties (darkness, fluorescence) with chemistry and weather data

Optical Microscopy


- Why do we use optical microscopy?
 - New approach at characterizing aerosol samples
 - Build upon other experimental work and correlate with XRF techniques and (later down the road...weather data)
- Nikon Fluorescence Microscope with CCD Camera

Optical Results

Ebam filters



July 18, 2004 10X



July 19, 2004 10X

Nikon Fluorescence Microscope with CCD Camera



Conclusion

- Titanium and Iron are present in aerosols in variable and sometimes high concentrations
- At first pass, optical darkness of filters correlates with weather characteristics suggesting that high humidity and rain events correlate with higher concentrations of metals

Further Work

- All data collected will be correlated to weather data, specifically humidity, precipitation and wind direction.
- We hope to be able to characterize aerosol chemistry by back-tracking to determine source.

References

- <http://www.atmosphere.mpg.de/enid/n2.html>